



White Paper

Andy Roebuck, Vitalograph


Defining the Lower Limits of Normal (LLN) in Spirometry

The use of percent predicted when assessing lung function is widely used and an 80% 'cut-off' for a lower limit of normal (LLN) is commonly utilised and considered during diagnosis and monitoring of respiratory disease. However, the 80% of predicted 'cut-off' for what may be considered 'normal' is not based on any scientific foundation and has been considered as statistically invalid³.

As many authors in the past have shown^{1,2,4}, scatter plots of data of respiratory parameters such as FEV1 show in healthy non-smokers, a proportional change of scatter as the value declines with age. As such, a LLN should follow the same pattern. When looking at the raw data scatter (Graph 1) using 80% predicted LLN shows a progressively increasing LLN as age increases with no statistical reasoning. Implementing this as a marker of abnormality, would create a large percentage of false-positive results in both young and elderly subjects. Additionally, with an ever increasing ageing population, this could be regarded as a potential problem when diagnosing and monitoring respiratory disease.

[Graph 1. The decline of a statistical 1.64SD LLN compared with an 80% predicted LLN.](#)

The comparison of the observed and predicted values which determines whether the difference




between the measurements is compatible with the scatter distribution about the predicted value, is strongly recommended. Such a method needs to:

- Determine how normal a result is.
- Measure how far the result is from the mean (predicted) value.
- Express it in terms of the number of standard deviations (SD) from that mean (figure 1)

Thus, enabling clearer statistical weighting and rationale during a time of no clinically proven categorisation of severity of disease as opposed to an arbitrary 80% LLN.

Figure 1. A bell-curve representing the LLN of 1.64 SD accepting a false positive of 5% and a lower false positive of 2.5% if LLN was increased to 1.96 SD.



- Figure 1 shows the prevalence of data that exceeded or fell below a predicted value (Z-score of 0) from the chosen regression equation.
- Both the ERS and ATS have endorsed the LLN as 1.645 SD below the mean value applying the 5th percentile LLN.
- The 'Z-score' or Standard Residual (SR) show the number of standard deviations a result is from the mean value.
- By implementing a recognised LLN of 1.64 SD, observations lower than -1.64 would occur in only 5% of the reflected population.
- Implementing a 1.96 SD LLN would imply that this was only observed in 2.5% of the 'healthy' reference population.
- Using this as a statistical marker for abnormality of a result in the absence of any other scientifically proven method, could have a beneficial effect in assisting with respiratory diagnosis and monitoring as well as showing strong application in research clinical trials for patient outcomes and recruitment when boundaries or 'cut-off' markers are used.

Table 1 shows that subject (a) had a measured FEV1 which would be considered 'normal' when implementing an 80% LLN, but shows that this result would be considered 'abnormal' as it is in the 5th percentile of their population normality. Conversely, subject (b) shows an FEV1 of 75% predicted but is within a statistical LLN (90th percentile) based on their given population normality.

Table 1.

Height (cm)	Age (Yrs)	FEV1 Obs (L)	FEV1 Pred (L)	FEV1 80% Pred (L)	Z-score
(a) 185	39	3.48	4.36	3.48	-1.72
(b) 177	70	2.34*	3.12	2.49	-1.52

Nb. *Subject (b) achieved 75% predicted FEV1 but still fell within a statistical LLN Z-score (not less than -1.64).

It is worth serious consideration the recommendation to implement the use of Z-score (SR) when considering diagnosis, monitoring and treating respiratory conditions. A non-proven LLN of 80% predicted for what is considered a 'normal' result may prove to be costly for both patient and practice.

References:

1. Quanjer PhH, Tammeling GJ, Cotes JE, Pedersen OF, Peslin R, Yernault JC. Lung volumes and forced ventilatory flows. *Eur Respir J* 1993; 6 suppl. 16: 5-40. Erratum *Eur Respir J* 1995; 8: 1629
2. Miller MR, Quanjer PH, Swanney MP, Ruppel G, Enright PL. Interpreting lung function data using 80% predicted and fixed thresholds misclassifies more than 20% of patients. *Chest* 2011; 139: 52-59
3. Sobol BJ, Sobol PG. Percent of predicted as the limit of normal in pulmonary function testing: a statistically valid approach. *Thorax* 1979; 34: 1-3
4. Quanjer PH, Stanojevic S, Cole TJet al.and the ERS Global Lung Function Initiative. Multi-ethnic reference values for spirometry for the 3-95 years age range: the Global Lung Function 2012 equations. *Eur Respir J* 2012; 40: 1324–1343