Peak Flow Meter

Universal For Adults & Children | Universal Für Erwachsene & Kinder | Universel pour adultes et enfants | Universale per adulti e bambini

CE **Model 4300**  Vitalograph, Gort Road Business Park, Ennis, Co. Clare, V95 HFT4, Ireland
EN ISO 23747
2797 Medical Device | Medizinisches Gerät | Dispositif médical | Dispositivo medico

GB Instructions for Use

Setting Action Zones

Zone markers may be set with a stylus by the healthcare professional to give personal boundaries for implementing the asthma action plan.

To Use Peak Flow Meter

1. Slide pointer to the bottom of the scale.
2. Hold device horizontally in front of your mouth.
3. Inhale as deeply as possible.
4. Seal your lips firmly around the mouthpiece.
5. Blow out as HARD and as FAST as possible for a second or longer.
6. Take care that tongue or teeth do not obstruct the airflow. A 'spitting' action or coughing will give false readings.
7. Take note of the reading.
8. Repeat steps 1 to 7 three times and record the highest reading.
9. Take action according to your personal asthma management plan.

Cleaning

Clean the outer surfaces after use. We recommend using an ordinary alcohol wipe (IPA 70%). Pay special attention to the mouthpiece area.

Indications For Use

Device measures peak expiratory flow (PEF) for use on adults and paediatrics, 5 years and older, in primary care, clinics, hospitals and at home, under the supervision of a healthcare professional.

Use In Clinic

When used for multiple subjects Vitalograph® intends that a mouthpiece adaptor (available as an accessory) and disposable Bacteria/Viral Filter (BVF™) be used. The Peak Flow Meter should be externally cleaned between each patient.
Cat. No 28501: EcoBVF – Bacterial/Vital Filter (100)

Disposal

The peak flow meter and its accessories constitute minimally soiled waste from human healthcare and should be disposed of in line with local requirements. Peak Flow Meter  ABS

Contraindications, Warnings, Precautions and Adverse Reactions


Note: Please read all this information before using this Peak Flow Meter.

- When used for multiple subjects Vitalograph intends that a new BVF be used for every subject to prevent cross contamination. Using a new BVF provides a significant level of protection of the subject, the device and the user against cross contamination. A BVF is for single use only.
- This device is intended to help you assess how well your asthma is under control, by measuring your peak flow which can then be monitored over time.
- Never attempt to dismantle the unit. This can cause faulty Peak Expiratory Flow scores.
- If used at home symptoms must take precedence over device measurements.
- If you think that the device is not reading correctly, advise your healthcare professional immediately.
- If the device is used for longer than it's specified life, the accuracy of the device may deteriorate.
- Store in a clean dry place.
- Any serious incident that has occurred in relation to the device should be reported to Vitalograph or its Authorized Representative and the Regulatory Authorities of the country.

Warranty

Your device is guaranteed for one year*. Replace if it is faulty, otherwise replace the unit every three years. *Excepting accidental/transit damage or inappropriate use of the device.

CE Notice

Marking by the  symbol indicates compliance of this device to the Medical Devices Directive of the European Community. Such marking is indicative that the device meets or exceeds the referenced technical standards.

FDA Notice

Caution: Federal Law restricts this device to sale by, or on the order of a physician.

Declaration of Conformity

Product: Vitalograph Model 4300

Vitalograph hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

European Medical Devices Directive 93/42/EEC, as amended.

This device is classified as 1m per Annex V per Annex VII section 5 of the MDD
Canadian Medical Device Regulation
{CMDR SOR/98-282}

FDA Quality System Regulation {QSR} 21 CFR 820;

EN ISO 13485. Medical devices. Quality management systems. Requirements for regulatory purposes

Certifying Body {For 93/42/EEC and CMDR}: British Standards Institute, Certificate Nos.

CE 00772, CE 85553, MD 82182.

Technical Specifications

Product: Peak Flow Meter

Model Number: 4300

Material: ABS plastic

Accuracy: ± 10L/min or ± 10% of the reading
(EN ISO 23747:2015)

Repeatability: ± 10L/min or ± 5% of the reading
(EN ISO 23747:2015)

Frequency: Profile A/B difference less than 15 L/min or 12% (EN ISO 23747:2015)

Highest resistance to flow: 0.00377 kPa/L/min @ 720 L/min

Measurement Range: 50-800 L/min BTPS

Operating Temperature Range: ATS/ERS: Limits: 17-37°C. **Design Limits:** 10-40°C

Peak Flow Meter performance standards:

EN ISO 23747:2015; ATS/ERS:2019

Designed & Manufactured Under: EN ISO 13485,

FDA 21CFR820, CMDR & JPAL

Size: 138.0mm x 56mm x 23mm

Weight: 78.4g (Net)

Package Contents: Peak Flow Meter, Instructions for Use

Frank Keane

Signed on behalf of Vitalograph (Ireland) Ltd.
Frank Keane, CEO

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Signé au nom de Vitalograph (Ireland) Ltd.