

VITALOGRAPH ALPHA SERIES

The Flexible All Round Spirometers For Use With or Without a PC



A Global Leader in Respiratory Diagnostics

Vitalograph ALPHA

Highly portable lightweight all-in-one desktop spirometer.
Easy to use. No PC required.

The Vitalograph® ALPHA™ is the spirometer of choice for respiratory care professionals within general healthcare and occupational health around the world due to its ease of use, portability, integral printer and long term reliability.

- Integral thermal printer
- Durable measurement technology without the need for expensive disposable sensors, turbines, spirottes or flow tubes
- Optimal hygiene via Bacterial Viral Filters (BVF™)
- Quick and easy calibration verification routine as recommended by good practice & international spirometry guidelines
- Simple and intuitive navigation
- Built-in interpretation of measurement results
- Mains or battery powered
- Includes Vitalograph Reports Software with the ability to create editable PDFs
- Features latest GLI spirometry predicted equations and Z-scores



Vitalograph ALPHA Touch

The Flexible All Round Spirometer For Use With or Without a PC

Vitalograph ALPHA Touch is an easy to use, lightweight desktop spirometer that offers the benefits of portability and a built-in printer. A large, high resolution colour touchscreen display, combined with an intuitive icon driven user interface, makes the ALPHA Touch quick and simple to operate.

- The ability to use with or without a PC
- Touchscreen for modern and intuitive navigation
- Patient management database for storing more than 10,000 patients
- Intuitive icon driven user interface
- Choice of child friendly incentive animations to help improve patient cooperation
- Optional upgrade to Spirotrac® Software enables automatic data exchange, via HL7 or GDT transfer protocols, with electronic medical record (EMR) systems



✓ Low Running Costs ✓ Optimal Hygiene

Ordering Information

6000 ALPHA
6000 ALPHA Touch
6000 ALPHA Touch (With 7000 Spirotrac Software)

Consumables and Optional Extras

28554 2820 Eco BVF with Bite Lip (75)
28553 Eco BVF with Bite Lip + Disposable Nose Clip (75)
36020 Precision Syringe - 3 Litre
66149 Thermal Printer Paper (5 rolls)

Vitalograph®, Vitalograph®, Spirotrac®, ALPHA™ and BVF™ are trademarks or registered trademarks of Vitalograph Ltd.

CE
2797 56026_3

**our services
cover the globe**

Vitalograph Inc.

Toll Free: 800 255 6626
Tel: (913) 730 3200
Fax: (913) 730 3232
contact@vitalograph.com
13310 West 99th Street
Lexen, Kansas, 66215,
USA

Authorised Dealer Stamp

Vitalograph Ltd, Hong Kong/China

sales@vitalograph.cn
P.O. Box 812
Shatin Central Post Office
Hong Kong

Vitalograph Ltd, International

Tel: +44 1280 827120
Fax: +44 1280 823302
sales@vitalograph.co.uk
Maids Moreton, Buckingham
MK18 1SW,
England

Vitalograph GmbH,

Tel: +49 40 547391-0
Fax: +49 40 547391-40
info@vitalograph.de
Rellinger Straße 64a
D-20257 Hamburg
Germany

Vitalograph Ltd, UK

Tel: 01280 827110
Fax: 01280 823302
sales@vitalograph.co.uk
Maids Moreton, Buckingham
MK18 1SW,
England

Vitalograph (Ireland) Ltd

Tel: +353 65 686 4100
Fax: +353 65 682 9289
sales@vitalograph.ie
Gort Road Business Park
Ennis, Co Clare, V95 HFT4,
Ireland

Vitalograph®

vitalograph.com