

3-L Precision Syringe

For peace of mind regarding the accuracy of your spirometer

The 3-L Precision Syringe from Vitalograph® allows quick, easy accuracy checks for all spirometers, as recommended by the ATS/ERS Standardization of Spirometry 2019 Update and required by the Spirometry standard ISO 26782:2009.

The 3-L Precision Syringe provides reliable verification of accuracy and can be used with both flow or volume sensor systems and is compatible with any manufacturers' equipment. An in-line filter adapter is included to allow use of Bacterial-Viral Filters for protection against contamination.

Technical Specifications

Model Number:	2040
Volume Measurement:	3L
Measuring Accuracy:	0.5%
Operating Temperature:	15–40°C
QA/GMP Standards:	ISO 13485, FDA 21CFR820
Gross Size:	500mm x 170mm x 125 mm
Gross Weight:	2.5kg
Storage Relative Humidity:	10%–95%



Vitalograph®, Vitalograph® is a registered trademark of Vitalograph Ltd.

CE 062987_1

Vitalograph®

A Global Leader in Respiratory Diagnostics